**1. Indication**  
**GR-17 temporary** is a resin for the generative production of temporary anterior dental restorations based on image projection systems ( $\leq 405$  nm).  
**GR-17.1 temporary It** is a resin for the generative production of long term temporary dental restorations (US: up to 30 days) based on image projection systems ( $\leq 405$  nm).

**2. Contraindication**  
**GR-17 temporary** and **GR-17.1 temporary It** restorations are contraindicated ...

1. ... if a patient is known to be allergic to any of the ingredients.
2. ... for bridges with more than one pontic.
3. ... for every application that is not part of the indication (see above).

**3. Requirements**

**Software – Information available from:**  
 DENTAL WINGS INC · 160 Rue St-Viateur E Suite 710, Montréal QC H2T 1A8 CANADA  
 3SHAPE A/S · Holmens Kanal 7, 1060 Copenhagen, Denmark

**Hardware (3D Printing) – Information available from:**  
 ASIGA · 2/19-21 Bourke Road, Alexandria, NSW 2015, Australia  
 Structo Pte Ltd · 60 Mac Pherson Road, #06-08A, The Siemens Centre, Singapore (348615)  
 MIIICRAFT · 27 Queen St. East Suite 1401 Toronto, Ontario M5C 2M6, Canada  
 DAZZ3D · Rm K-2F, Building A, No. 9 Zhongshanyuan Rd., Shenzhen, China 518052

**Hardware (Post Curing) – Information available from:**  
 pro3dure medical (CD-1, CD-2) · Am Burgberg 13, D-58642 Iserlohn, Germany  
 Structo Pte Ltd · 60 Mac Pherson Road, #06-08A, The Siemens Centre, Singapore (348615)  
 NK Optik (Otoflash G171) · Isarstraße 2, 82065 Baierbrunn, Germany  
 Prusa Research a.s. (Prusa CW1) · Partyzánská 188/7a, 170 00 Praha 7, Czech Republic

**4. Material**  
**GR-17 temporary** and **GR-17.1 temporary It** consists of functional methacrylic resins and inorganic fillers with particle sizes from 0.4 to 3 microns.

**5. Geometric presetting**  
**Connector areas (minimum):**  
 anterior bridges 12 mm<sup>2</sup>  
 posterior bridges 14 mm<sup>2</sup>

**Minimum wall thickness:**  
 occlusal 1.5 mm (central fissure)  
 circumferential 1 mm

**6. Material parameters**  
 Radiation penetration depth controllable by exposure time  
 $\geq 100 \mu\text{m}$  (4 mils)  
 $\geq 125 \mu\text{m}$  (5 mils)

**7. Manufacturing process (fig. 1–10, page 1)**

1. Prepare data (CAD & build preparation).
2. Choose process parameters (build style, etc.).
3. Transfer prepared data to 3D printing device.
4. Prepare 3D printing – shake bottle.
5. Fill resin tank of 3D printing device.
6. Build the parts.
7. Clean parts (with IPA  $\geq 97\%$  or equivalent cleaning agent) approx. 4 min. in an ultrasonic bath or equivalent device – precleaning recommended.
8. Dry parts (with compressed air until no residues of IPA or equivalent cleaning agent).
9. Post curing (10 min.): inert atmosphere recommended (use adequate light curing device).
10. Finish parts.

**8. Individualization and cementation**  
 The temporary restorations can be individualized by light curing composite materials. The reduction of the restoration should not exceed 0.3 mm (incisal and vestibular area). Please follow the instructions of the light curing composite manufacturers. For cementation common provisional dental cements (free of Eugenol) and adhesives are suitable. The inner areas of the crown should be conditioned by sandblasting or other commonly used techniques. Please follow the instructions and indications of the corresponding manufacturers.

**9. Notice**  
 Please follow the instructions provided by the software manufacturer in terms of parameter settings and design recommendations. Please follow the instructions provided by the hardware manufacturer in terms of parameter settings/printing & post curing recommendations. To avoid detrimental effects on material quality do not expose the liquid material to irradiation under any circumstances. Deviations from the described manufacturing processes or storing conditions may lead to different mechanical and optical properties of the material. Ensure personal protective gear during processing.  
**Caution:** Polymerized resins are chemically resistant – avoid stains on clothing! Avoid any contact with skin and eyes! In case of accidental contact, rinse with adequate running water, consult a doctor if necessary. The lot number and the best before date are indicated on each packaging. In case of claims please always indicate the lot number of the product. Do not use the product after expiry of the best before date. Dispose of contents/container according to official regulations.

**10. Hazard statements**  
 Causes skin irritation. May cause an allergic skin reaction. Causes serious eye damage. May cause respiratory irritation. Harmful to aquatic life with long lasting effects.

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**Distributor (US):**  
 pro3dure medical LLC  
 9825 Valley View Road · Eden Prairie, MN 55344  
 Phone: 952-426-1928 · Fax: 952-681-7515  
 info@pro3dure.com · www.pro3dure.com · Made in Germany

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