Upcoming entrained bubbles get eliminated through a holding time of approx. 1 hour with FotoDent® cast.

When embedding FotoDent® cast, please note the manufacturer's data of the used embedding materials. To ensure that the resin combusts without any residue, a holding time of 60 minutes at 800–900 °C is recommended.

To avoid damage of material quality do not expose the liquid material to light under any circumstances. Variations from the described manufacturing process may lead to modified mechanical characteristics and colour variations of the FotoDent® cast material.

Possible dangers

Please note the safety instructions for processing FotoDent® cast.

Hazard statements: **H315** Causes skin irritation. **H319** Causes serious eye irritation. **H317** May cause an allergic skin reaction. **H335** May cause respiratory irritation.

Precautionary statements: **P261** Avoid breathing dust/fume/gas/mist/vapours/spray. **P280** Wear protective gloves/protective clothing/eye protection/face protection. **P304+P340** IF INHALED: Remove person to fresh air and keep comfortable for breathing. **P305+P351+P338** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. **P405** Store locked up. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Contains: Diurethandimethacrylate, triethylene glycol dimethacrylate, phenyl bis(2,4,6-trimethylbenzoyl)-phosphineoxide

Lot number / Expiry date

The lot number and the expiry date are indicated on each FotoDent® cast packaging. In case of claims please always indicate the lot number of the product. Do not use the product after the expiry date.

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Intended use

FotoDent® cast is a light curing resin for the production of cauterisation profiles for casting technology using 385 nm / 405 nm LED based stereolithography. For post curing of FotoDent® cast built parts, the light curing units PCU LED N₂ and FotoDent® flash are recommended, which are optimal balanced to the production process.

Processing

- During the application of FotoDent® cast we recommend wearing personal protective equipment e. g. suitable gloves, safety goggles etc.
- FotoDent® cast should be well shaken approx. 1 hour before use.
- Carefully pour FotoDent® cast into the pre-determined container of the production unit.
- Remove all bubbles with a cleaned object.
- Processing temperature:

Rapid Shape D and DII printer	23 ± 3 °C
Asiga MAX UV	35 ± 3 °C
- Select the product specific parameter settings for your process.
- When the building process is finished a direct post treatment is recommended.
- After the platform is taken up a dripping off time of approx. 10 minutes is recommended.
- Resin coated parts get clean with isopropanol (97 %) within a time of about 10 minutes. For a particular effective cleaning with isopropanol the use of explosion-protected ultrasonic units is recommended.
- Post curing **PCU LED N₂**:

Rapid Shape process	6 min 80 % light power (vacuum)
Asiga process	6 min 80 % light power (vacuum)
- Post curing **FotoDent® flash**:

Rapid Shape process	8 min+ 1 min gas pre-flow
Asiga process	10 min + 1 min gas pre-flow
- Impurity due to the construction or a break of the material and impurity due to operation mistakes cannot be excluded. Thanks to the low viscosity however it is possible to filtrate FotoDent® cast. It is recommended to regularly take out the container of the production unit, to homogenise and filtrate the content.

Dreve Dentamid GmbH is not liable for any damages caused by improper application of the material.

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Dreve Dentamid GmbH haftet nicht für Schäden, die durch fehlerhafte Anwendung des Materials hervorgerufen werden.