

sagemax.com

sagemax®





Sintering

NexxZr® T, NexxZr®+

	Number of units	Duration	Phase	Temperature	Heating/ Cooling rate	Temperature	Heating/ Cooling rate	Holding time
					°C/min.		°F/min.	
Standard	1-5	~ 3.7	1	20-1300	30	68-2372	54	30
			2	1300-1530	40	2372-2786	72	60
			3	1530-900	15	2786-1652	27	-
			4	900-80	20	1652-176	36	-
	5-10	~ 5.2	1	20-1300	30	68-2372	54	60
			2	1300-1530	40	2372-2786	72	120
			3	1530-900	15	2786-1652	27	-
			4	900-80	20	1652-176	36	-
Long	1-20	10.7	1	20-900	10	68-1652	18	10
			2	900-1530	3	1652-2786	5.4	150
			3	1530-80	8	2786-176	14.4	-
	>21	11.7	1	20-900	10	68-1652	18	10
			2	900-1530	3	1652-2786	5.4	210
			3	1530-80	8	2786-176	14.4	-
Over night	unlimited	~ 14.3	1	20-250	2	68-482	3.6	-
			2	250-1530	4	482-2786	7.2	240
			3	1530-80	8	2786-176	14.4	-



sagemax[®]

Sintering

NexxZr[®] + Multi, NexxZr[®] T Multi

	Number of units	Duration	Phase	Temperature	Heating/ Cooling rate	Temperature	Heating/ Cooling rate	Holding time
				°C	°C/min.	°F	°F/min.	min.
Speed	1-5	4.8	1	20-1000	60	68-1832	108	10
			2	1000-1530	3	1832-2786	5.4	60
			3	1530-1100	50	2786-2012	90	-
			4	1100-100	60	2012-212	108	-
Long	unlimited	9.6	1	20-900	10	68-1652	18	30
			2	900-1500	3.3	1652-2732	5.94	120
			3	1500-900	10	2732-1652	18	-
			4	900-300	8.3	1652-572	14.94	-

 Sagemax Bioceramics, Inc., 34210 9th Ave. South, Suite 118
Federal Way, WA 98003, USA

 See instructions **CE 0123** Rx only. For dental use only!



sagemax.com



sagemax®

Zirconia Portfolio

NexxZr®



Product	NexxZr ⁺ Multi	NexxZr ⁺	NexxZr T Multi	NexxZr T
Strength	Incisal 650 MPa Dentin 850 MPa	Preshaded 880 MPa White 1000 MPa	Incisal 650 MPa Dentin 1200 MPa	1270 MPa
Translucency	Incisal 49% Dentin 45%	46%	Incisal 49% Dentin 41%	42%
Colors	A1, A2, A3, B1, B2, C2, D2	White, A1, A2, A3, A3.5, B1, B2	BL1, BL2, BL3, BL4, A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4	White, BL1, BL2, BL3, A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4
Geometries	W-98, Z-95	W-98, Z-95	W-98, Z-95	W-98, Z-95
Thicknesses	16, 20 mm	10, 12, 14, 16, 18, 22, 25 mm	16, 20, 25 mm	10, 12, 14, 16, 18, 22, 25 mm

sagemax.com

S01049/EN/2023-05/Rev. 1